



3D Ball Mill (3D Reactor)
PC control model
(Program & Manual operation possible)

May, 2018




NAGAO SYSTEM

Chassis color has been changed from gray to white. Gray color is also possible.



PC control model (Program & Manual operation possible)



Speed changeable in increments of 10 rpm.



Pattern operation is also possible.
Variable time, hold time and rotation speed can
be adjusted.

回転パターン設定 Version 1.0.0.1

パターン選択 --- 新規追加 ---

パターン名 有効 読込

STEP1~5 STEP6~10 繰り返し回数
0

	モーター 1			モーター 2		
	変動時間(秒)	保持時間(秒)	回転速度(rpm)	変動時間(秒)	保持時間(秒)	回転速度(rpm)
STEP 1	0	0	0	0	0	0
STEP 2	0	0	0	0	0	0
STEP 3	0	0	0	0	0	0
STEP 4	0	0	0	0	0	0
STEP 5	0	0	0	0	0	0

登録 更新 削除 閉じる

It is also possible to compare measured rpm with program mode.



Rotation pattern can be set.

